# **Trophen**

Tropicamide 0.8% & Phenylephrine HCI 5% Sterile Eye Drops

#### **Composition:**

Trophen Sterile Eye Drops: Each ml contains Tropicamide USP 8 mg and Phenylephrine

Hydrochloride USP 50 mg

Preservative: Benzalkonium Chloride Solution BP 0.02%.

Vehicle: HPMC USP 0.5%.

#### Indications:

As a mydriatic, Trophen is indicated in some pre & postoperative states and for several ophthalmologic examinations like during ophthalmoscopy, slit-lamp examination, retinal photography, laser treatment, or adjunct in the treatment of anterior uveitis. Also it may be used in temporary lowering of intraocular pressure in glaucoma.

#### Dosage and administration:

Instill 1 drop in the eye(s) 15-20 minutes before examination. If examination is not conducted within 20-30 minutes, an additional drop may be placed in the eye(s) to prolong the effect.

#### **Contraindications:**

Known hypersensitivity to any ingredient of the preparation, and narrow angle glaucoma.

## Side effects:

Ocular side effects include transient stinging and raised intra-ocular pressure; on prolonged administration, local irritation hyperemia, edema and conjunctivitis may occur. Systemic effects include arrhythmias, hypertension, and coronary artery spasm.

# **Drug Interaction:**

Phenylephrine may interact with systemically administered Monoamino Oxidase Inhibitors (MAOIs).

## **Precaution:**

Children and elderly; cardiovascular diseases, tachycardia, hypertension, diabetes.

# Storage:

- Store in a cool & dry place, protect from light.
- Do not refrigerate.
- Do not use for longer than one month after first opening the bottle.
- Keep out of the reach of children.

#### Packing:

Trophen Sterile Eye Drops: Each Plastic dropper bottle contains 5 ml/10 ml sterile eye drops.



